

## PRODUCT INFORMATION

Orange Serum Agar

Cat. No. O15-101

### DESCRIPTION

Orange Serum Agar is used for the cultivation of aciduric microorganisms, particularly those associated with the spoilage of citrus products. The low pH of fruit juices makes citrus fruit products susceptible to spoilage by yeasts, molds, *Lactobacillus* and *Leuconostoc*.

### FORMULA (g/L)

Orange Serum Extract	13.5 g	Tryptone	10.0 g
Yeast Extract	3.0 g	Dextrose	4.0 g
Dipotassium Phosphate	2.5 g	Agar	17.0 g

Final pH: 5.5 ± 0.2 at 25 °C

\*Grams per liter may be adjusted or formula supplemented to obtain desired performance.

### PREPARATION

Mix 25 grams of the medium in one liter of purified water until evenly dispersed. Heat with repeated stirring and boil for one minute to dissolve completely. Distribute and autoclave at 121°C for 15 minutes.

### QUALITY CONTROL SPECIFICATIONS

1. The powder is homogenous, free flowing and tan to beige.
2. Visually the prepared medium is trace to slightly hazy and light to yellowish-beige.
3. Expected cultural response after 2-5 days at 30 °C.

ORGANISM	RESULT
<i>Aspergillus brasiliensis</i> ATCC 16404	Good Growth
<i>Lactobacillus fermentum</i> ATCC 9338	Good Growth
<i>Leuconostoc mesenteroides</i> ATCC 23386	Good Growth
<i>Saccharomyces cerevisiae</i> ATCC 9763	Good Growth



## **STORAGE**

Store the sealed bottle containing the dehydrated medium at 2 to 30°C. Once opened and recapped, place the container in a low humidity environment at the same storage temperature. Protect it from moisture and light. The dehydrated medium should be discarded if it is not free flowing or if the color has changed from the original color.

Version 01 – Date 10/04/23