



PRODUCT INFORMATION

Peptic Digest of Animal Tissue

Cat. No. P16-102

DESCRIPTION

Peptic Digest of Animal Tissue, otherwise known as Meat Peptone, is used for the cultivation of microorganisms. It provides nitrogen, amino acids and carbon for microbiological culture media.

TYPICAL PHYSICAL & CHEMICAL CHARACTERISTICS

TEST TYPE	SPECIFICACION
Appearance	Beige powder
Amino Nitrogen	≥ 3.4%
Total Nitrogen	≥ 10.0%
Ash	< 15.0%
pH (2% solution)	6.5 – 7.5

STORAGE

Store the sealed bottle containing the dehydrated medium at 2 to 30°C. Once opened and recapped, place the container in a low humidity environment at the same storage temperature. Protect it from moisture and light. The dehydrated medium should be discarded if it is not free flowing or if the color has changed from the original color.

Version 01 – Date 09/08/23